



Open Interface

DASware® migrate software—the online matchmaker in bioprocess control

DASware® software solutions are well-known as a synonym for advanced bioprocess development.



Add DASware® to Your Existing Bioprocess Controller

DASware migrate provides the gateway to use all DASware solutions and DASGIP® modules with legacy benchtop bioreactor control units. Users of Eppendorf New Brunswick™ and third-party controllers from Sartorius®, Applikon®, and others will profit from adapting the excellent DASware software functionalities to their bioprocesses. Process development according to QbD guidelines, automation by integration of online analytics, and interconnectivity with corporate IT systems results in time- and cost-effective time-to-market.


Users benefit from a comfortable, user-friendly and comprehensive bioprocess management tool for multiple bioreactors on the market.



DASware migrate has been proven effective in the integration of Eppendorf New Brunswick BioFlo®/ CelliGen® 115 & 310, BioFlo 320 and 415, and CelliGen® BLU as well as Sartorius BIOSTAT® Bplus DCU and Applikon 1010, 1030 & ez-Control. It also allows for mixed system control with user-friendly batch-to-batch comparison.

Process Development in Line with QbD Guidelines

DASware migrate facilitates Quality by Design (QbD) approaches in bioprocess development, regardless which legacy benchtop bioreactor equipment is used. It opens up the possibility to use all functions of DASware software:



DASware® control: Advanced bioprocess control

Our bioprocess control software DASware control follows our easy »Point-Click-Grow« concept. Readily applicable recipes with incorporated batch functionality are designed to reflect process control needs. Stored as recipe templates, these predefined strategies can be emailed, loaded, changed, saved, and executed immediately.

DASware® access:

Remote control—anytime from anywhere

Running bioreactor processes can be remotely observed and controlled from everywhere. All modifications to process values are tracked and logged in a central database for complete documentation whereas data security is ensured by established IT security mechanisms. DASware access uses Wi-Fi, Intranet, VPN, and 3G connections to provide access via PC and Notebook. The universal DASGIP iApp supports access from iPhone®, iPod touch® and iPad®.

DASware® design: Supports Design of Experiments

With DASware design automated DoE workflows can be applied to the bioprocess to identify critical process parameters (CPPs) systematically. An integrated full factorial DoE builder as well as integration of a large variety of third-party DoE designs e.g. JMP® and multivariate analysis tools facilitates finding a proper design space for manufacturing processes. A recipe generator supports multiple system set-ups. Following the Point-Click-Grow concept they can be carried out on a set of bioreactors with a single mouse click.

DASware® connect: Integration with legacy control

Using DASware connect all bioprocess information can be integrated with legacy control systems and corporate historians such as Emerson® DeltaV, Siemens® SIMATIC PCS 7, ABB® 800 xA, OSIsoft® PI System and Matrikon® OPC Historian facilitating a company-wide access to all relevant bioprocess data. Also, it supports interfacing with scientific software packages like LabVIEW® and MATLAB®.

DASware® analyze:

Integration of third-party analyzers and autosamplers

An OPC network protocol allows for interconnectivity between bioreactors systems, autosamplers, and analyzers such as cell counters, biomass monitors, nutrient analyzers, mass spectrometers or HPLCs. DASware analyze makes direct feedback from the bioreactor system possible in response to online measured analytical data. It supports online calculations as well as event- and data-driven decisions.

DASware® discover:

Comprehensive information management

Configurable and retrievable CPPs can be online or retrospectively added to process runs. Intuitive queries allow for near real-time retrieval of runtime information from a SQL Server® database. A Microsoft® Excel® report generator provides recipe information, process information, and event reporting. Users can simultaneously compare process information from either current or historical runs.



DASGIP® modules in use with legacy equipment

DASware migrate easily supports third-party bioreactor controllers with DASGIP GA4 exhaust analyzers, DASGIP OD4 biomass monitors, DASGIP MP8 multi pumps and DASGIP MX4/4 massflow-controlled gas mixing stations.

> Move your bioprocess to the next level—with DASware software solutions.
www.ependorf.com/DASware

Ordering information

| Description | Order no. |
|---|-------------|
| DASware® migrate , license for 1 vessel | |
| for New Brunswick™ controllers | 76DWMIGNB |
| for third-party systems | 76DWMIGTP |
| DASware® migrate , incl. PC, OS | |
| package to operate New Brunswick™ systems (up to 16 vessels) | 76DWMIGNBPC |
| package to operate third-party systems (up to 16 vessels) | 76DWMIGTPPC |
| DASware® control , incl. PC, OS, and licenses | |
| for 4-fold DASGIP® system | 76DGCS4 |
| for 8-fold DASGIP® system | 76DGCS8 |
| for 16-fold DASGIP® system | 76DGCS16 |
| for 4-fold DASbox® system | 76DXCS4 |
| for 8-fold DASbox® system | 76DXCS8 |
| for 16-fold DASbox® system | 76DXCS16 |
| for 24-fold DASbox® system | 76DXCS24 |
| DASware® access , remote access support (web and iApp) for 1 vessel | 76DWACC |
| DASware® analyze , OPC client standard (OPC DA e.g. for ext. analyzer), for 1 vessel | 76DWANA |
| DASware® analyze , license for serial/Modbus® integration (e.g. for ext. biomass sensors), for 1 vessel | 76DWANAM |
| DASware® analyze , OPC client professional incl. 1x tunneller lic. (OPCDA e.g. for ext. analyzer with autosampler) | |
| for 4 vessels | 76DWANA4P |
| for 8 vessels | 76DWANA8P |
| for 12 vessels | 76DWANA12P |
| DASware® analyze , cable and license | |
| for 4 Aber Futura® sensors | 76DWANA4AF |
| for 4 Hamilton® Fogale sensors | 76DWANA4HF |
| DASware® connect , OPC server (OPC DA for ext. PCS), for 1 vessel | 76DWCON |
| DASware® design , DoE and local information management, license for 1 vessel | 76DWDOE |
| DASware® discover Client License , for 1 vessel (SQL Server®-based information management) | 76DWDIS |
| DASware® discover Information Management Server , PC hardware, OS software, and server licence | 76DWDISPC |
| DASware® discover Server License , SQL Server®-based information management | 76DWDISS |

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · 22331 Hamburg · Germany
 eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Sartorius® and BIOSTAT® are registered trademarks of Sartorius AG, Germany. Applikon® is a registered trademark of Applikon B.V., The Netherlands. Microsoft®, Excel®, and SQL Server® are registered trademarks of Microsoft Corporation, USA. LabVIEW® is a registered trademark of National Instruments Corporation, USA. MATLAB® is a registered trademark of The Mathworks, Inc., USA. Emerson® is a registered trademark of Emerson Electric Co., USA. Siemens® is a registered trademark of Siemens AG, Germany. ABB® is a registered trademark of ABB Asea Brown Boveri Ltd., Switzerland. OSIsoft® is a registered trademark of OSIsoft, LLC, USA. Matrikon® is a registered trademark of Matrikon Inc., Canada. JMP® is a registered trademark of SAS institute Inc, USA. iPhone®, iPad®, and iPod touch® are registered trademarks of Apple Inc., USA. Modbus® is a registered trademark of Schneider Electric USA, Inc. Eppendorf® and the Eppendorf logo are registered trademarks of Eppendorf AG, Germany. New Brunswick™ is a trademark of Eppendorf AG, Germany. DASGIP® and DASware® are registered trademarks of DASGIP Information and Process Technology GmbH, Germany. CelliGen® and BioFlo® are registered trademarks of Eppendorf, Inc., USA.

All rights reserved, including graphics and images. Copyright ©2016 by Eppendorf AG.
 Order no. A765 111 020/GB2/1T/0116/EBC/STEFF · Carbon neutrally printed in Germany.